
blue●dot®

Pulse Oximeter

Model: B-PO252 BT

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
Contents

Included in the Package	2
Brief Description of Front Panel	2
General Description	5
Measurement Principle	6
Safety Information	8
Components	19
Product Features	19
Intended Use	20
Setup	22
Operation Instructions	23
Battery Installation	25
Cleaning and Disinfection	26
Maintenance and Storage	28
Specification	29
Graphical Plot of Data Points	33
Assessing the Validity of a SpO ₂ Reading	35
Troubleshooting	37
Appendix 1 EMC Information	41
Possible Problems and Resolutions	51
Symbol Definitions	54


Included in the Package

1. Pulse Oximeter * 1
2. Pouch Bag * 1
3. AAA battery * 2
4. Manual * 1
5. Lanyard * 1



Brief Description of Front Panel

The pulse plethysmogram displays corresponding with the user's pulse beat. The height of the plethysmogram shows the user's pulse strength.

The oximeter supports multi-directional display, and the screen can display measured value, pulse wave, and bar graph with 4 directions display.

(Signal inadequacy display interface)

When the oxygen saturation and pulse rate signals are incomplete, the measured values are not displayed.

Caution:

The bar graph acts as an indicator of signal incompleteness, the higher the bar graph the stronger the pulse. When the bar graph

is not beating smoothly, the accuracy of the measured value may be reduced at this time. When the bar graph tends to be beating smoothly, the measured value read is the best value.

General Description

SpO₂ stands for peripheral capillary oxygen saturation. Oxygen saturation is defined as the ratio of oxyhemoglobin (HbO₂) to the total concentration of hemoglobin (i.e. Oxyhemoglobin + reduced hemoglobin) present in the blood. It is an important physiological parameter involved in respiration and circulation. The Pulse Oximeter feature herein is small, portable, non-invasive and easy to use. The user only needs to insert a finger into the chamber to measure his/her SpO₂ and pulse rate. The oximeter is suitable for home healthcare environments.

Declaration


EMC of this product complies with IEC60601-1-2 standard. The materials which the user can come into contact with have no toxicity and no action on tissues comply with ISO10993-1, ISO10993-5 and ISO10993-10.

Measurement Principle

Oxygenated blood absorbs light preferentially at 905 nm (near infrared light), whereas deoxygenated blood absorbs light preferentially at 660nm (red light). A pulse oximeter works by passing a beam of red and infrared light through a pulsating capillary bed and then measure the amount of red and infrared light emerging from the tissues via a sensor. To improve accuracy, the Oximeter uses a proprietary algorithm to collect data from

pulsatile arterial blood and excludes local noise from the tissues. The relative absorption of light by oxyhemoglobin (HbO) and deoxyhemoglobin is then calculated according to the Beer-Lambert's law and a quantitative measurement of the users' oxyhemoglobin status i.e. oxygen saturation level (SpO_2) is derived.

Due to the sensitivity of the pulse oximeter, finger should be kept stationary during measurement. It is recommended that you use this device for measurement before or after sports. Do not use for continuous monitoring.

Safety Information

DANGER

Indicates an imminent hazard that, if not avoided, will result in death or serious injury.

Dangers

There are no dangers that refer to the product in general.

Cautions

1. This device does not have alarm function.
2. Prolonged use or user conditions may require periodic changes to the sensor position. Change sensor position at least every 2 hours and check skin integrity, circulation condition, and correct alignment.

3. Although the part expected to come into contact with the human body has been biologically evaluated and the biological safety meets the standard requirements, very few people may have allergic reactions, and patients with allergic reactions should stop using it immediately.
4. The accuracy of pulse oximetry measurements may be adversely affected by high ambient light. Please avoid direct light (such as surgical lights or sunlight) during measurements.
5. Do not disassemble, install, or repair the host without authorization. Do not look directly at the light-emitting devices, which may cause eye damage.
6. The maximum temperature of the blood oxygen probe and the contact surface of the human body does not exceed 41°C, which is measured by the temperature measuring instrument.

7. The product's blood oxygen waveform has been normalized.
8. The data update period is less than 30s, and the data averaging and signal processing affects the display of SpO₂ and pulse rate: The blood oxygen (SpO₂) and pulse rate (PR) values displayed on the oximeter are the values obtained by calculating and averaging multiple sets of data acquired within a set period of time. The data is calculated once for each pulse wave, and the data is a single real-time blood oxygen and pulse rate value. The data averaging method is related to the pulse rate value; when the pulse rate value is below 50 bpm, the last 12 seconds of data are averaged, when the pulse rate value is between 50 bpm and 120 bpm, the last 8 seconds of data are averaged, and when the pulse rate value is greater than 120 bpm, the last 6 seconds of data are averaged. The displayed blood oxygen and pulse rate values are displayed in 1 second update intervals.

9. The product is not recommended for use during exercise and is not suitable for continuous monitoring of patients.
10. Not used for hypoperfusion oximetry. Veins beat rhythmically.

WARNING

Indicates a potential hazard or unsafe practice that, if not avoided, will result in death or serious injury.

Warnings

1. Federal (U.S.) law restricts this device to sale by or on the order of a physician.
2. Explosive hazard — DO NOT use the equipment in the environment with tinder such as anesthetic.
3. DO NOT use the equipment while the patient is being scanned by MRI or CT.
4. Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.
5. Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in

improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

6. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
7. The equipment is NOT intended for neonate and infant, and the patient's finger thickness should be between 8 to 25.4 mm.
8. It is recommended that the equipment should be inspected before use, when there is obvious damage, stop using the equipment.

9. Discomfort or pain may appear if using the equipment ceaselessly, especially for microcirculation barrier patients, it is recommended that the equipment should not be used on the same finger more than 10 minutes.
10. For some patients who need a more careful examination of the measurement site, the oximeter shall not be placed in edema or fragile organization.
11. The equipment is just a clinical diagnosis auxiliary equipment. The physiological data displayed on the equipment are for reference only and cannot be directly used for diagnostic interpretation.
12. It is not recommended to use the equipment in high frequency environment such as electrosurgical equipment.

13. Do not have the equipment immersed in liquid.
14. Prevent children from swallowing the equipment or its accessories. Children must be accompanied by adult guardian when using products.
15. Please follow local ordinances and recycling instructions regarding disposal or recycling of the equipment and batteries.
16. No maintenance or servicing the product when using.
17. Users are NOT permitted to maintain the equipment by themselves.
18. There are NO replaceable components in the equipment.

CAUTION

Indicates a potential hazard or unsafe practice that, if not avoided, could result in minor personal injury or product/property damage.

Cautions

The equipment is designed to measure the percentage of arterial oxygen saturation of functional hemoglobin. Factors that may degrade pulse oximeter performance or affect the accuracy of the measurement include the following:

- Do not apply the pulse oximeter on the same arm as a blood pressure cuff, arterial catheter or infusion line(s) (IVs)
- Excessive light, such as sunlight or direct home lighting

- Moisture in the equipment
- Finger is beyond recommended size range
- Poor pulse quality
- Venous pulsations
- Anemia or low hemoglobin concentrations
- Cardio green and other intravascular dyes
- Carboxyhemoglobin
- Methemoglobin
- Dysfunctional hemoglobin
- Artificial nails or fingernail polish
- Excessive movement during testing may result in incorrect readings
- Improper wearing position may result in incorrect reading

The light (the infrared is invisible) emitted from the equipment is harmful to the eyes, so the user and the maintenance man should not look at the light.

NOTE

1. The patient is an intended operator, and all the functions can be safely used.
2. Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.

Components

The oximeter consists of probe, electronic circuits, and display and plastic enclosures.

NOTES

- The probe is the hole in the middle of the equipment to which the finger insert.
- The probe is the Applied Part of the equipment.

Product Features

- Lightweight, portable and easy to use.
- OLED screen shows pulse rate, SpO₂, perfusion index.
- Large font display
- Low battery indicator
- Auto shut down if no signal is detected within 15 seconds.

Intended Use

The pulse oximeter is a reusable device and intended for spot-checking of oxygen saturation and pulse rate for use with the finger of adult patients.

Intended User

Pulse Oximeter could be used by Healthcare Professional and lay person.

Intended Patient

Pulse Oximeter is used to adult patients.

Indications

The pulse oximeter can effectively estimate and predict the blood oxygen concentration in the periphery of the body by non-invasive operation, so as to achieve the role of monitoring the changes of the respiratory disease with using the finger pulse oxygen as the monitoring index.

Contraindications

It is not intended to be used when the following happens:
pulse oximeter should not be fixed on the tissue injury site;
Not used for blood oxygen monitoring in hyperkinetic patients;
Not used for blood oxygen measurement under hypoperfusion.

Set up

After the oximeter is turned on, press the power on button for a second to enter the setting interface.

Press the power on button once to switch among the settings, keep pressing the button for a second to change the current settings of sounds on or off and other parameters.

Remind Setup	*
Sound Reminder	on
Beep	off
Demo	off
Restore	ok
Brightness	4
Exit	

Limit Setup	*
SpO2 Hi	100
SpO2 Lo	94
PR Hi	130
PR Lo	50
+/-	+
Exit	

1. A sound is sounded when the blood oxygen level falls below the set threshold to alert the user

2. A sound is sounded when the pulse rate exceeds the set threshold to alert the user

Note: The sound function is only used to remind the user to check their blood oxygen level and is not used for any alarm purposes. Any time users have doubts about their blood oxygen saturation, they should check their blood oxygen measurement in time.

Operation Instructions

1. Install two AAA batteries into battery compartment correctly.
2. Insert one of your fingers into the finger chamber of the pulse oximeter.

Note: The fingernail should be facing the chamber which contains the sensor. Finger should also be inserted completely into the chamber. Otherwise, measurement will be inaccurate.

3. Press the power-on button to turn on the pulse oximeter.
4. Finger and body should not tremble during measuring.
5. Read correct data from displayed screen.

Notes:

- When your finger is inserted into the oximeter, your nail surface must be on the sensor.
- The results may be wrong if you do not place your finger thoroughly in the oximeter.
- Please use medical alcohol to clean the silicon which touches the finger inside of oximeter, and use alcohol to clean the test finger before and after each test. (The silicon inside of the oximeter belongs to medical silicon, which has no toxins and does no harm to the skin).

Battery Installation

Put the two AAA batteries into battery compartment in correct polarities.

Notes:

- Battery polarities should be correctly installed. Otherwise, damage may be caused to the device.
- Please remove the batteries if the Oximeter will not be used for a long time.
- This product should be kept away from children, because small parts such as battery cases can cause children to choke.

Cleaning and Disinfection

1. Cleaning

Your equipment should be cleaned on a regular basis. If there is heavy pollution or lots of dust and sand in your place, the equipment should be cleaned more frequently. Recommended cleaning agents are:

- Ethanol (70%)

To clean your equipment, follow these rules:

- Shut down the pulse oximeter.
- Clean the display screen using a soft, clean cloth dampened with a glass cleaner.
- Clean the exterior surface of the equipment and probe using a soft cloth dampened with the cleaner.
- Wipe off all the cleaning solution with a dry cloth after cleaning if necessary.
- Dry your equipment in a ventilated and cool place.

To avoid damage to the equipment, follow these rules:

CAUTIONS

- Always dilute according the manufacturer's instructions or use lowest possible concentration.
- Do not immerse the equipment in liquid.
- Do not pour liquid onto the equipment or accessories.
- Never use abrasive materials (such as steel wool or silver polish), or erosive cleaners (such as acetone or acetone-based cleaners).
- If you spill liquid onto the equipment, contact us or your service personnel.

2. Disinfection

Clean the pulse oximeter before disinfecting it.

The recommended disinfectant is ethanol 70%. Disinfection steps are the same as cleaning. Do not disinfect the instrument by using high-temperature/high-pressure disinfecting gas.

Never use ETO or formaldehyde for disinfection.

Maintenance and Storage

1. Replace the batteries in time when low battery indicator is lighted.
2. Clean surface of the fingertip oximeter before it is used in diagnosis for users.
3. Remove the batteries inside the battery cassette if the oximeter will not be operated for a long time.
4. It is best to preserve the product in a place where ambient temperatures is $-20^{\circ}\text{C} \sim 50^{\circ}\text{C}$ and relative humidity is 10% - 95%.
5. It is recommended that the product should be kept in a dry environment anytime. A wet ambient might affect its lifetime and even might damage the product.

6. Avoid exposure to direct sunlight.
7. Avoid excessive radioactive infrared rays or ultraviolet ray
8. Please follow the law of the local government to deal with_S.
used battery.

Specification

1. Display Type: OLED
2. SpO₂:
 - Measurement range: 0% - 100%
 - Display range: 35% - 100%
 - Resolution: 1%
 - Accuracy: $\pm 2\%$ for 70% - 100%; less than 70% is unspecified.

3. Pulse Rate:
Measurement range: 25 bpm - 250 bpm
Resolution: 1 bpm
Accuracy: ± 2 bpm
Pulse Intensity: Bargraph Indicator
4. Perfusion Index:
Measurement range: 0% - 50%
Resolution: 1%
5. Power Requirements: Two AAA alkaline Batteries
Power consumption: 30 mA (Normal)
Low battery indication: When the battery voltage is $1.9V \pm 0.2V$, the oximeter will automatically shut down
Battery Life: Two AAA 1.5V, 600 mAh alkaline batteries could be continuously operated as long as 24 hours.
6. Dimension: Length x Width x Height:
63.4 mm x 36.0 mm x 33.6 mm
Weight: ~38 g (without battery)

7. Environment Requirements:

Temperature: 10°C - 40°C

Humidity (non-condensing): 15% - 95%

Atmospheric pressure : 70 kPa - 106 kPa

8. Storage/transportation conditions

Temperature: -20°C - 50°C

Humidity (non-condensing): 10% - 95%

Atmospheric pressure: 70 kPa - 106 kPa

9. Measurement Performance in Low Perfusion Condition: 0.3%

10. Safety Classification

Type of protection against electric shock:

Internally powered equipment

Degree of protection against electric shock:

Type BF-Applied part (non-defibrillation proof)

Operating mode: Continuous operation

Degree of protection against hazards of explosion: IP22

11. Range of the peak wavelengths: RED: 660 nm/IR: 905 nm
12. Maximum optical output power: 1.2 mW
13. Use-life: Two years
14. Communication mode: Bluetooth 4.2

Note:

1. SpO₂ and pulse rate accuracy are the essential requirements.
2. There is no alarm that includes the capability to detect an SpO₂ or pulse rate physiological alarm condition is available.
3. Data Update duration is not applicable since it is not intended to be connected with other device to display the reading.
4. The oximeter is for spot checking, frequency of use is not applicable.

Graphical Plots of Data Points


Figure 1 Bland-Altman graph for SpO₂-SaO₂


Figure 2 The regression line of the measured data

Assessing the Validity of a SpO₂ Reading

You can check the quality of the pleth wave and the stability of the SpO₂ values on the Patient Monitor to assess whether the sensor functions properly and whether the SpO₂ readings are valid. Always use these two indications simultaneously to assess the validity of a SpO₂ reading.

NOTE:

1. The SpO₂ accuracy has been validated in human studies against arterial blood sample reference measured with a CO-oximeter. Pulse oximeter measurements are statistically distributed, only about two-thirds of the measurements can be expected to fall within the specified accuracy compared to CO-oximeter measurements. The volunteer population in the studies composed of local healthy men and women from age

18 to 46, with variations of skin pigmentations. And the accuracy cannot be assessed by a function tester.

2. The pulse rate accuracy is obtained by comparison to the pulse rate generated with an arterial oxygen simulator (also an electronic pulse simulator).
3. Generally, the quality of the SpO_2 pleth wave reflects the quality of the light signals obtained by the sensor. A wave of poor quality manifests a decline of the signal validity. On the other hand, the stability of the SpO_2 values also reflects the signal quality. Different from varying SpO_2 readings caused by physiological factors, unstable SpO_2 readings are resulted from the sensor's receiving signals with interference. The problems mentioned above may be caused by patient movement, wrong sensor placement or sensor malfunction. To obtain valid SpO_2 readings, try to limit patient movement.

Troubleshooting

Trouble	Possible Reason	Solution
The equipment can't be turned on.	The battery is drained away or almost drained away.	Please replace battery.
	The battery installation is incorrect.	Install the battery over again.
	The malfunction of the equipment.	Please contact the local service center.

The display is off suddenly.	The equipment is set to shut down automatically in 15 seconds when there is no correct physiological signals.	Normal
	The battery is almost drained away.	Please replace battery.
The SpO ₂ and Pulse Rate are not displayed stably.	The thoroughfare from photo detector to light emitting diode was sheltered by some objects.	Check and clean the inside of the probe especially the two windows of sensors.
	The finger is shaking or the patient is moving.	The patient tries to keep still.

	The finger is not placed inside deep enough.	Place the finger properly and try again.
	The finger's size is too big or too small.	Select the correct size finger to measure.
	Excessive ambient light	Avoid the excessive ambient light irradiation.
	Pulse rate value of the cyclical fluctuations.	The measurement is normal, and the patient has arrhythmia.

The SpO ₂ and Pulse Rate can't be displayed normally.	The finger is not properly positioned.	Place the finger properly and try again.
	The patient's SpO ₂ is too low to be detected.	Try again, go to a hospital for a diagnosis if you are sure the equipment works all right.
Pulse sound can't be turned off.	The key is bad.	Check the key and press again.
	The press time is not right.	Make sure the press time is 2~3 seconds.

Appendix 1 EMC Information

Guidance and manufacturer's declaration- Electromagnetic emission		
Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Model B-PO252 BT should assure that it is used in such an environment.		
Emissions	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Pulse Oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

RF emissions CISPR 11	Class B	The Pulse Oximeter is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2	N.A.	N.A.
Voltage fluctuations/flicker emissions IEC61000-3-3	N.A.	

Guidance and manufacturer's declaration - Electromagnetic immunity

The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of Pulse Oximeter should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.

Electrical fast transient/burst IEC 61000-4-4	Not applicable	Not applicable	
Surge IEC 61000-4-5	Not applicable	Not applicable	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	Not applicable	

Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC61000-4-6	Not applicable	Not applicable	
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	
NOTE: UT is the a.c. mains voltage prior to application of the test level			

Guidance and manufacturer's declaration - electromagnetic Immunity

The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of Pulse Oximeter should assure that it is used in such an environment.

Radiated RF IEC 61000-4-3 (Test specifications for ENCL OSU RE POR	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
	385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1,8	0.3	27	27

T IMM UNIT Y to RF wirel ess com muni cation s equi pme nt)	450	430 – 470	GMRS 460, FRS 460	FM ± 5 kHz deviati on 1 kHz sine	2	0.3	28	28
	710	704 – 787	LTE Band 13, 17	Pulse modula tion 217 Hz	0,2	0.3	9	9
	745							
	780							
	810	800 – 960	GSM 800/90 0, TETRA 800, iDEN	Pulse modula tion 18 Hz	2	0.3	28	28
	870							

	930		820, CDMA 850, LTE Band 5					
	1720	1700- 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modula tion 217 Hz	2	0.3	28	28
	1845							
	1970							

	2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modula tion 217 Hz	2	0.3	28	28
	5240	5100 – 5800	WLAN 802.11 a/n	Pulse modula tion 217 Hz	0,2	0.3	9	9
	5500							
	5785							

Guidance and manufacturer's declaration - electromagnetic Immunity

Radiated RF IEC61000-4-39 (Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields)	Test Frequency	Modulation	IEC 60601-1-2 Test Level (A/m)	Compliance level (A/m)
	30 kHz	CW	8	8
	134,2 kHz	Pulse modulation 2.1 kHz	65	65
	13,56 kHz	Pulse modulation 50 kHz	7,5	7,5

Possible Problems and Resolutions

Problems	Possible reason	Solution
The Oximeter fails to display the blood oxygen saturation levels and/or pulse rate.	1. The finger is not placed between the sensor and the Light Emitting Diode. 2. The user's blood perfusion is too low to be detected.	1. Make sure that the finger is placed right in between the sensor and the Light Emitting Diode. 2. Make sure nothing is restricting your blood flow.

SpO₂ or PR is shown unstably	1. The finger is not placed between the sensor and the Light Emitting Diode. 2. The user is moving his/her finger and/or body.	1. Make sure that the finger is placed right in between the sensor and the Light Emitting Diode. 2. Try to stay still during measurement.
The Oximeter can not be powered on	1. The batteries are drained. 2. The batteries are Incorrectly installed. 3. The Oximeter is defective and/or damaged.	1. Replace the batteries. 2. Please refer to the section 'Battery Installation' on how to install the battery correctly. 3. Contact the product owner.

The screen are suddenly off	1. If the Oximeter does not detect any signal within 15 sec, it will automatically power off. 2. The batteries are drained.	1. This is normal. Just turn on the pulse oximeter again. 2. Replace the batteries.
------------------------------------	--	--

There are no user-serviceable parts inside the oximeter. The cover should only be removed by qualified service personnel. If you are uncertain about the accuracy of any measurement, check the user's vital signs by alternate means; then make sure the oximeter is functioning correctly.

Do not spray, pour, or spill any liquid on the oximeter, its accessories, connectors, switches, or openings in the enclosure as this may damage the oximeter.

Symbol Definitions

Symbol	Definition
	Type BF Applied part
	Caution
% SpO2	Oxygen saturation
 /Min	Pulse rate
	The equipment should be sent to the special agencies according to local regulation for separate collection after its useful life.
	Alarm inhibit
	Refer to instructions manual

IP22	2 Protected against solid foreign objects of 12.5 mm Ø and greater ; 2 Protected against vertically falling water drops when ENCLOSURE tilted up to 15°
	Date of Manufacture
RoHs	RoHS mark
	CE mark
	Manufacturer
	Authorized representative in the European
	Medical device

	General symbol for recovery
LOT	Batch code

Note: The illustration used in this manual may differ slightly from the appearance of the actual product.

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Product Owner: Lifebox Asia Company Limited